



# Wisconsin Department of Justice Toxicology Unit Procedure Manual TXPM 6.5 Benzodiazepines/Z-drugs

## Section 6 – Analytical Methods, Specific

### Part 5 – Benzodiazepines/Z-drugs

#### Purpose and Scope

This document presents methods to analyze specimens suspected of containing benzodiazepines, z-drugs and various metabolites. These drugs are extracted from biological samples using solid phase extraction (SPE) technology and then analyzed by GC/MS or LC/MS/MS. If samples are analyzed by GC/MS, derivatization is required for some analytes. Urine samples are hydrolyzed prior to extraction and analysis. A specialized extraction procedure for clonazepam is also included.

#### Table of Contents

<b>Purpose and Scope</b> .....	1
<b>6.5.1 Definitions</b> .....	1
<b>6.5.2 Chemicals and Reagents</b> .....	1
<b>6.5.2.1 Chemicals</b> .....	1
<b>6.5.2.2 Reagents</b> .....	2
<b>6.5.3 Equipment and Supplies</b> .....	2
<b>6.5.4 Procedure</b> .....	3
<b>6.5.4.1 Specimen Criteria</b> .....	3
<b>6.5.4.2 Certified Reference Materials</b> .....	3
<b>6.5.4.3 Calibrators</b> .....	3
<b>6.5.4.4 Controls</b> .....	4
<b>6.5.4.5 Hydrolysis Procedures</b> .....	4
<b>6.5.4.5.1 Hydrolysis Procedure – Option 1- GC/MS</b> .....	4
<b>6.5.4.5.2 Hydrolysis Procedure – Option 2 – LC/MS/MS</b> .....	5
<b>6.5.4.6 Benzodiazepine Procedure Option 1 - SPE with TMS Derivatization and GC/MS Analysis</b> .....	5
<b>6.5.4.7 Clonazepam Procedure – SPE with TBDMS Derivatization</b> .....	7
<b>6.5.4.8 Benzodiazepine Procedure Option 2 – SPE with LC/MS/MS Analysis</b> .....	7
<b>6.5.4.9 LC/MS/MS Performance Characteristics</b> .....	9
<b>6.5.4.10 LC/MS/MS Instrument Parameters</b> .....	10
<b>6.5.4.11 Limitations of LC/MS/MS Analysis</b> .....	15
<b>6.5.4.12 LC/MS/MS Stability</b> .....	15
<b>6.5.5 Results, Conclusions, and Data Analysis</b> .....	16
<b>6.5.6 Reporting</b> .....	16
<b>6.5.7 References</b> .....	16
<b>6.5.8 Safety</b> .....	16
<b>6.5.9 Method Validation</b> .....	17

#### 6.5.1 Definitions

**Hydrolysis** – a process in which the conjugate bond is broken to allow the freed moiety to be detected.

**Total drug concentration** – The sum of the unbound/free and bound/conjugated drug.

#### 6.5.2 Chemicals and Reagents

##### 6.5.2.1 Chemicals

Abalonase ultra

Acetonitrile

BGTurbo enzyme (β-Glucuronidase)

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	1 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*



# Wisconsin Department of Justice Toxicology Unit Procedure Manual TXPM 6.5 Benzodiazepines/Z-drugs

BSTFA+  
Ethyl acetate  
 $\beta$ -Glucuronidase solution (IMCSzyme 3S)  
 $\beta$ -Glucuronidase G7017 (from *Helix pomatia*)  
 $\beta$ -Glucuronidase solution, from abalone  
MTBSTFA+  
Hexanes  
Instant buffer I  
Rapid hydrolysis buffer  
Methanol  
Methanol (LC/MS grade)  
Acetonitrile (LC/MS grade)  
Isopropanol (LC/MS grade)  
Methylene chloride  
Formic acid (LC/MS grade)  
LC/MS/MS water

## 6.5.2.2 Reagents

Abalonase enzyme solution  
0.1M Acetate buffer, pH 5.0  
0.1M Acetate buffer, pH 4.5  
1M Acetic acid, pH 2.4  
1 M Acetic acid: acetonitrile  
Alprazolam extraction solvent  
EA with 2% Ammonium Hydroxide  
Basic elution solvent  
Benzo extraction solvent  
Borate buffer, saturated, pH 9.0  
Clonazepam extraction solvent  
Enzyme hydrolysis reagent (from limpets or from *Helix pomatia*)  
0.1M Phosphate buffer, pH 6.0  
0.1M Phosphate buffer, pH 6.0 (LC/MS/MS)  
LC/MS/MS mobile phase - Methanol with 0.1% formic acid  
LC/MS/MS mobile phase - Water with 0.1% formic acid  
LC/MS/MS mobile phase - Water with 2% formic acid  
LC/MS/MS needle wash  
LC/MS/MS seal wash

See [TXPM 2.1 Chemicals and Reagents](#) for necessary chemicals and preparation instructions.

## 6.5.3 Equipment and Supplies

Centrifuge, capable of  $\geq 4,000$  rpm  
Evaporator/Concentrator  
GC/MS Instrument, (Agilent or similar) with applicable software, a compatible computer, and printer  
Heating Block with NIST traceable thermometer  
LC/MS/MS Instrument, (Waters Inc. or similar) with applicable software, a compatible computer, and printer  
PTFE septa and caps, pre-slit and in-house assembled  
Pipets – Pasteur, disposable  
Pipets – Transfer, disposable  
Pipettors – air displacement with disposable tips

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	2 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*



# Wisconsin Department of Justice

## Toxicology Unit Procedure Manual

### TXPM 6.5 Benzodiazepines/Z-drugs

- Pipettors – positive displacement (QP specifications) with disposable tips
- Repeater pipette
- Screw Caps – Teflon® lined
- Solid Phase Extraction Manifold
- SPE Columns, Clean Screen® Extraction Column (ZSDAU020) from UCT (200 mg)
- Test Tube Rocker or Rotator
- Tubes – round/conical bottom borosilicate glass with screw tops or disposable culture tubes
- Vials for LC/MS/MS Autosampler with inserts and pre-slit caps
- Vials for GC Autosampler with inserts and caps
- Vortex Mixer

#### 6.5.4 Procedure

##### 6.5.4.1 Specimen Criteria

- 1 mL whole blood, biological fluids, or tissue homogenates with preservatives, if applicable, for GC/MS analysis.
- 500 µL whole blood or urine with preservatives, if applicable for LC/MS/MS analysis.
- Samples will be removed from storage and brought to room temperature before analysis.

##### 6.5.4.2 Certified Reference Materials

Spiking solutions will be prepared from certified reference materials. Spiking solutions for qualitative analytes may be prepared from reference materials (see [TXPM 2.4 Reference Materials and Certified Reference Materials](#)).

Internal standard solution may be deuterated versions of the analytes.

##### 6.5.4.3 Calibrators

Calibrators are prepared utilizing the spiking solutions. See Table 1 for an example of calibrator concentrations for GC/MS analysis. See Table 2 for an example of calibrator concentrations for LC/MS/MS analysis.

Table 1 – An example of calibrator concentrations for GC/MS.

Drug or Metabolite	Cal A ug/L	Cal B ug/L	Cal C ug/L	Cal D ug/L
Diazepam	500	250	100	50
7-Aminoflunitrazepam	100	50	20	10
Midazolam	100	50	20	10
Estazolam	100	50	20	10
Alprazolam	100	50	20	10
Triazolam	100	50	20	10
7-Aminoclonazepam	100	50	20	10
α-OH-alprazolam	100	50	20	10
α-OH-triazolam	100	50	20	10
Nordiazepam	500	250	100	50
Oxazepam	200	100	50	25
Lorazepam	100	50	20	10
Temazepam	100	50	20	10
2-OH-ethylflurazepam	100	50	20	10

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	3 of 17



# Wisconsin Department of Justice Toxicology Unit Procedure Manual TXPM 6.5 Benzodiazepines/Z-drugs

Table 2 – An example of calibrator concentrations for LC/MS/MS.

Drug or Metabolite	Cal 7 ug/L	Cal 6 ug/L	Cal 5 ug/L	Cal 4 ug/L	Cal 3 ug/L	Cal 2 ug/L	Cal 1 ug/L
7-aminoclonazepam	240	160	80	40	20	10	5
Alprazolam	240	160	80	40	20	10	5
Clonazepam	240	160	80	40	20	10	5
Diazepam	1000	500	200	100	40	20	10
Etizolam	240	160	80	40	20	10	5
Lorazepam	240	160	80	40	20	10	5
Nordiazepam	1000	500	200	100	40	20	10
Oxazepam	1000	500	200	100	40	20	10
Temazepam	1000	500	200	100	40	20	10
Zolpidem	240	160	80	40	20	10	5
$\alpha$ -OH-alprazolam	240	160	80	40	20	10	5

#### 6.5.4.4 Controls

Negative and positive controls must meet requirements defined in [TXPM 4.4 Quantitative and Qualitative Quality Controls](#) or [TXPM 4.9 Liquid Chromatographic and Tandem Mass Spectral Quality Control](#).

#### 6.5.4.5 Hydrolysis Procedures

Urine specimens should be hydrolyzed to determine the total drug concentration prior to extraction using one of the following hydrolysis procedures. Blood specimens may also be hydrolyzed for GC/MS. A glucuronide positive control shall be used to ensure that the hydrolysis was effective, if available.

##### 6.5.4.5.1 Hydrolysis Procedure – Option 1- GC/MS

- 1) Into appropriately labeled tubes add:
  - Hydrolysis reagent (per Table 3).
  - Spike internal standard and controls, as necessary.
  - 1 mL unknown or blank matrix.
- 2) Gently mix the tubes.
  - For bloods, allow to equilibrate for 10 minutes prior to hydrolysis.
- 3) Cap and if applicable, incubate (per Table 3).
- 4) If applicable, cool the tubes.
- 5) To each tube add:
  - 7 mL 0.1M acetate buffer, pH 5.0
- 6) Continue with step 7 in TXPM 6.5.4.6 below.

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	4 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*



**Wisconsin Department of Justice  
Toxicology Unit Procedure Manual  
TXPM 6.5 Benzodiazepines/Z-drugs**

Table 3 – Hydrolysis Methods

Hydrolysis Method	Hydrolysis Reagent	Minimum Time	Temperature (+/- 5°C)
Enzyme hydrolysis reagent (from Helix pomatia or limpets)	1mL enzyme hydrolysis reagent	2 hours	60°C
Campbell Science β-glucuronidase solution, from abalone	≥60μL β-glucuronidase, from abalone 200μL acetate buffer, pH 5.0	30 minutes	65°C
Kura Biotech BGTurbo High Efficiency Recombinant β-Glucuronidase	400μL Instant Buffer I 200μL BGTurbo Enzyme 1100 μL H <sub>2</sub> O	No wait	Room temperature
United Chemical Technologies (UCT) Abalonase Ultra β-glucuronidase enzyme	1mL Abalonase Enzyme Solution	30 minutes	65°C
Integrated Micro-Chromatography Systems (IMCS) 3S genetically modified β-glucuronidase	400μL Rapid Hydrolysis Buffer 20μL β-glucuronidase	30 minutes	55°C

**6.5.4.5.2 Hydrolysis Procedure – Option 2 – LC/MS/MS**

- 1) Into appropriately labeled tubes add:
  - 200 μL Instant Buffer I
  - 100 μL BGTurbo
  - 550 μL H<sub>2</sub>O
  - Spike internal standard and controls, as necessary
  - 500 μL unknown or blank matrix
- 2) Gently mix the tubes.
- 3) Incubate at room temperature for zero minutes (i.e. no wait hydrolysis)
- 4) Add 2 mL 0.1M phosphate buffer, pH 6.0 (LC/MS/MS).
- 5) Continue with step 5 in TXPM 6.5.4.8 below.

**6.5.4.6 Benzodiazepine Procedure Option 1 - SPE with TMS Derivatization and GC/MS Analysis**

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	5 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*



# Wisconsin Department of Justice

## Toxicology Unit Procedure Manual

### TXPM 6.5 Benzodiazepines/Z-drugs

- 1) Label tubes accordingly.
- 2) Spike internal standard, calibrators and controls as necessary.
- 3) May add 200-500 uL H<sub>2</sub>O to each tube.
- 4) Pipet 1 mL blank matrix into calibrators and controls and 1 mL unknown into appropriately labeled tubes.
- 5) Vortex for 5-10 seconds and equilibrate for  $\geq 10$  minutes.
- 6) To each tube add:
  - 8 mL 0.1M acetate buffer, pH 5.0
- 7) Cap and rock for  $\geq 10$  minutes.
- 8) Centrifuge for  $\geq 10$  minutes.
- 9) Set up the SPE manifold with SPE columns, Clean Screen® Extraction Column (ZSDAU020).
- 10) Prep the SPE columns with:
  - 1 mL methylene chloride
  - 3 mL methanol
  - 2 mL H<sub>2</sub>O
  - 1 mL 0.1M acetate buffer, pH 5.0
- 11) Transfer the supernatant from the tubes into the appropriate SPE column. Avoid transferring any of the pellet from the bottom of the tube.
- 12) Sequentially wash the columns with:
  - 1 mL H<sub>2</sub>O
  - 3 mL 1M acetic acid, pH 2.4
- 13) Dry the columns high pressure for  $\geq 12$  minutes.
- 14) Turn off the high pressure and wash the columns with:
  - 1 mL hexanes
- 15) Re-dry the SPE columns under high pressure for  $\geq 5$  minutes.
- 16) Wash the columns with:
  - 4 mL ethyl acetate
- 17) Briefly turn on the high pressure to draw out the last of the solvent from the columns.
- 18) Re-dry the SPE columns under high pressure for  $\geq 3$  minutes.
- 19) Turn off the high pressure. Put the labeled tubes into the manifold and elute with:
  - 3 mL basic elution solvent.
- 20) Briefly turn on the high pressure to draw out the last of the solvent out of the columns.
- 21) Evaporate to dryness at  $\leq 50^\circ\text{C}$  under nitrogen.
- 22) To each tube add:
  - 50 uL acetonitrile
- 23) Vortex and transfer to appropriately labeled autosampler vials with inserts. Reserve tubes for step 27. Cap the vials with in-house assembled/pre-slit **PTFE** septa caps.
- 24) The samples are now ready for underivatized GC/MS analysis.  
*Note: Underivatized benzodiazepines may include diazepam, 7-aminoflunitrazepam, alprazolam, estazolam, midazolam, and triazolam.*
- 25) After underivatized GC/MS analysis: to each autosampler vial perform a quantitative solvent rinse by doing the follow steps **3** times:
  - Add 100 uL acetonitrile
  - Vortex
  - Transfer to the appropriate tube
- 26) Evaporate to dryness at  $\leq 50^\circ\text{C}$  under nitrogen.
- 27) To each tube add:
  - Blood samples with 25 uL acetonitrile and 25 uL BSTFA+
  - Urine samples with 25 uL acetonitrile and 75 uL BSTFA+
- 28) Cap and incubate at  $90\pm 5^\circ\text{C}$  for  $\geq$  one hour.
- 29) Centrifuge for  $\geq 3$  minutes.
- 30) Transfer the solution into the appropriately labeled autosampler vials with inserts. Cap the vials with new in-house assembled/pre-slit **PTFE** septa caps.

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	6 of 17

This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.



# Wisconsin Department of Justice

## Toxicology Unit Procedure Manual

### TXPM 6.5 Benzodiazepines/Z-drugs

The samples are now ready for analysis by GC/MS.

*Note: Derivatized benzodiazepines may include 7-aminoclonazepam, 7-aminoflunitrazepam, alpha-hydroxyalprazolam, alpha-hydroxytriazolam, nordiazepam, oxazepam, lorazepam, temazepam and 2-hydroxyethylflurazepam.*

#### 6.5.4.7 Clonazepam Procedure – SPE with TBDMS Derivatization

- 1) Label tubes appropriately.
- 2) Spike internal standard, calibrators and controls as necessary.
- 3) May add 200-500  $\mu$ L H<sub>2</sub>O to each tube.
- 4) Pipet 1 mL blank matrix into calibrators and controls and 1 mL unknown into appropriately labeled tubes.
- 5) Vortex for 5-10 seconds and equilibrate for  $\geq$ 10 minutes.
- 6) To each tube add:
  - 8 mL 0.1M phosphate buffer, pH 6.0
- 7) Cap and rock for  $\geq$ 10 minutes.
- 8) Centrifuge for  $\geq$ 10 minutes.
- 9) Set up the SPE manifold with SPE columns, Clean Screen® Extraction Column (ZSDAU020).
- 10) Prep the SPE columns with:
  - 3 mL methanol
  - 2 mL H<sub>2</sub>O
  - 1 mL 0.1M phosphate buffer, pH 6.0
- 11) Transfer the supernatant from the tubes into the appropriate SPE column. Avoid transferring any of the pellet from the bottom of tube.
- 12) Wash the columns with:
  - 3mL clonazepam extraction solvent
- 13) Dry the columns under high pressure for  $\geq$ 12 minutes.
- 14) Wash the columns with:
  - 3 mL hexanes
- 15) Dry the columns under high pressure for  $\geq$ 12 minutes.
- 16) Turn off the high pressure. Put labeled tubes into the manifold and elute with:
  - 2 aliquots of 3 mL ethyl acetate.
- 17) Briefly turn on the high pressure to draw out the last of the solvent from the columns.
- 18) Evaporate to dryness at  $\leq$ 50°C under nitrogen.
- 19) To each tube add:
  - 25  $\mu$ L ethyl acetate
  - 25  $\mu$ L MTBSTFA+
- 20) Vortex, and transfer into appropriately labeled autosampler vials with inserts. Cap the vials.
- 21) Derivatize at room temperature for at least 20 minutes.

The samples are now ready for analysis by GC/MS.

#### 6.5.4.8 Benzodiazepine Procedure Option 2 – SPE with LC/MS/MS Analysis

- 1) Label tubes accordingly.
- 2) Spike internal standard, calibrators and controls as necessary.
- 3) Pipet 500  $\mu$ L blank matrix into calibrators and controls and 500  $\mu$ L unknown into appropriately labeled tubes.
- 4) Add 3 mL 0.1M phosphate buffer, pH 6.0 (LC/MS/MS) to each tube.
- 5) Vortex for 5-10 seconds and equilibrate for  $\geq$ 5 minutes.
- 6) Centrifuge for  $\geq$ 10 minutes.
- 7) Set up the SPE manifold with SPE columns, Clean Screen® Extraction Column (ZSDAU020).

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	7 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*



# Wisconsin Department of Justice

## Toxicology Unit Procedure Manual

### TXPM 6.5 Benzodiazepines/Z-drugs

- 8) Prepare the SPE columns with:
  - 3 mL methanol
  - 3 mL H<sub>2</sub>O (LC/MS/MS grade)
  - 1 mL 0.1M phosphate buffer, pH 6.0 (LC/MS/MS)
- 9) Transfer the supernatant from the tubes into the appropriate SPE column without transferring the pellet. Load at 0-1 mL/min ( $\leq 5$  psi on PPM), no faster to avoid analyte breakthrough.
- 10) Sequentially wash the columns with:
  - 3 mL H<sub>2</sub>O (LC/MS/MS grade)
    - Ensure no blood is on the walls of the column. Options include using a kimwipe, cotton tip applicator, or using the water in this step to rinse the walls.
  - 2 mL 1M acetic acid:acetonitrile
    - Ensure minimal to no aqueous is on the walls of the column
  - 2 mL hexanes
- 11) Dry the columns under high pressure for  $\geq 20$  minutes.
- 12) Put the labeled tubes into the manifold and elute with:
  - 3 mL EA with 2% Ammonium Hydroxide
    - Collect eluent with minimal airflow (1 drop per second approximately).
- 13) Evaporate to dryness at  $\leq 40^{\circ}\text{C}$  under nitrogen.
- 14) Return the labeled tubes into the manifold and elute with:
  - 3 mL EA with 2% Ammonium Hydroxide
    - Collect eluent with minimal airflow (1 drop per second approximately).
- 15) Briefly turn on the high pressure to draw out the last of the solvent out of the columns.
- 16) Evaporate to dryness at  $\leq 40^{\circ}\text{C}$  under nitrogen.
- 17) To each tube:
  - add 40  $\mu\text{L}$  methanol (LC/MS grade)
  - Vortex
  - Transfer the solution into appropriately labeled autosampler vials
- 18) Repeat step 17.
- 19) To each autosampler vial:
  - Add 120  $\mu\text{L}$  H<sub>2</sub>O (LC/MS/MS grade)
- 20) Cap the vials with pre-slit PTFE septa caps.

The vials are ready for LC/MS/MS analysis

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	8 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*





# Wisconsin Department of Justice Toxicology Unit Procedure Manual TXPM 6.5 Benzodiazepines/Z-drugs

## 6.5.4.9 LC/MS/MS Performance Characteristics

Table 4 - Quantitative analytes

Drug or Metabolite	LOD (µg/L)	LLOQ (µg/L)	Working range (µg/L)	Dilution* (1:1) 250 µL	Dilution* (1:4) 100 µL	Dilution* (1:9) 50 µL
7-Amino clonazepam	5	5	5-240	Yes	Yes	Yes
Alprazolam	5	5	5-240	Yes	Yes	Yes
a-OH alprazolam	5	5	5-240	Yes	Yes	Yes
Clonazepam	5	5	5-240	Yes	Yes	Yes
Diazepam	10	10	10-1000	Yes	Yes	Yes
Etizolam	5	5	5-240	No	No	No
Lorazepam	5	5	5-240	Yes	Yes	Yes
Nordiazepam	10	10	10-1000	Yes	Yes	Yes
Oxazepam	10	10	10-1000	Yes	Yes	Yes
Temazepam	10	10	10-1000	Yes	Yes	Yes
Zolpidem	5	5	5-240	Yes	Yes	Yes

\*When an unknown requires dilution to report quantitative results, a dilution QPC must be prepared the same as the unknown and run with the same dilution factor as the unknown.

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	9 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*



# Wisconsin Department of Justice Toxicology Unit Procedure Manual TXPM 6.5 Benzodiazepines/Z-drugs

Table 5 - Qualitative analytes  
*Limit of detection is 5 µg/L for the below analytes*

Drug/Metabolite	Drug/Metabolite
2-OH-ethylflurazepam	Estazolam
7-Aminoflunitrazepam	Flualprazolam
8-Aminoclonazepam	Flubromazepam
Adinazolam	Flubromazolam
a-OH-clonazolam	Flunitrazepam
a-OH-etizolam	Flurazepam
a-OH-flubromazolam	Lormetazepam
a-OH-triazolam	Midazolam
Bromazepam	Norchlordiazepoxide
Bromazolam	Nimetazepam
Chlordiazepoxide	Nitrazepam
Clobazam	Phenazepam
Clonazolam	Pyrazolam
Delorazepam	Triazolam
Demoxepam	Zaleplon
Diclazepam	Zopiclone

#### 6.5.4.10 LC/MS/MS Instrument Parameters

**Column:**

Waters Acquity UPLC BEH C18 1.7 µm 2.1 x 100 mm  
Optional: Pre-column or inline filter unit

**LC Parameters:**

Sample Manager Temperature Setpoint: 4°C ±12.0 °C  
Column Temperature Setpoint: 50 °C ±5.0 °C  
Injection Volume: 1 µL  
Flow rate: 0.400 mL/min

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	10 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*



# Wisconsin Department of Justice Toxicology Unit Procedure Manual TXPM 6.5 Benzodiazepines/Z-drugs

### Example mobile phases:

Mobile phase A: LC/MS/MS water

Mobile phase B: Methanol (LC/MS grade)

Mobile phase C: Water with 2% formic acid

Table 6 – Quaternary Solvent Manager Gradient Example

Time (minutes)	% A	%B	%C	Curve
Initial	65	30	5	6
2.50	30	65	5	6
3.25	25	70	5	6
4.50	18	77	5	6
4.51	5	90	5	6
4.80	5	90	5	6
5.91	65	30	5	6
9.00	65	30	5	6

Note: Variations in gradient may exist due to instrument capabilities, column properties, etc

### Example mobile phases:

Mobile phase A: Water with 0.1% formic acid

Mobile phase B: Methanol with 0.1% formic acid

Table 7 – Binary Solvent Manager Gradient Example

Time (minutes)	% A	%B	Curve
Initial	65	35	6
0.76	65	35	6
3.26	30	70	6
4.01	25	75	6
5.26	18	82	6
5.27	5	95	6
6.56	5	95	6
6.67	65	35	6
9.76	65	35	6

Note: Variations in gradient may exist due to instrument capabilities, column properties, etc

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	11 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*



**Wisconsin Department of Justice  
Toxicology Unit Procedure Manual  
TXPM 6.5 Benzodiazepines/Z-drugs**

**MS/MS Parameters:**

Capillary Voltage Setpoint: 0.7 kV  
 Ionization: ESI, positive  
 Source Temperature Setpoint: 150 °C  
 Desolvation temperature Setpoint: 500 °C  
 Cone Gas Flow Setpoint: 20 L/hr  
 Desolvation Gas Flow Setpoint: 800 L/hr

Note: Variations in capillary voltage may exist between different instruments

Table 8 - Ion Transitions (**Quantitation ions in bold**)

Compound	Precursor ion	Product Ions	Cone (V)	Collision (eV)
Nordiazepam	271	<b>140</b>	18	28
		165		
D5-Nordiazepam	276	<b>140</b>	59	30
		165		26
Nitrazepam	282	180	45	36
		207		34
7-Aminoflunitrazepam	284	135	2	28
		227		24
Diazepam	285	<b>154</b>	2	26
		193		32
7-Aminoclonazepam	286	<b>121</b>	14	30
		222		24
Norclordiazepoxide	286	227	34	22
		232		34
Oxazepam	287	104	48	40
		<b>241</b>		24
Demoxepam	287	104	52	22
		180		20
D4-7 Aminoclonazepam	290	<b>121</b>	28	32
		226		24
D5-Diazepam	290	154	4	28
		<b>198</b>		32
D7-7 Aminoflunitrazepam	291	230	92	30
		<b>138</b>	90	26
D5-Oxazepam	292	<b>246</b>	38	18
		274		14
Estazolam	295	205	6	38
		267		24
Nimetazepam	296	221	40	34
		268		22
Chlordiazepoxide	300	192	25	30
		227		24
D5-Estazolam	300	210	14	42
		272		24

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	12 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*



**Wisconsin Department of Justice  
Toxicology Unit Procedure Manual  
TXPM 6.5 Benzodiazepines/Z-drugs**

Temazepam	301	177	24	40
		<b>255</b>	14	16
Clobazam	301	224	40	34
		259		20
Delorazepam	305	140	59	30
		206		34
D5-Temazepam	306	<b>177</b>	38	40
		198		34
Zaleplon	306	236	50	26
		264		20
Zolpidem	308	<b>235</b>	58	40
		263		25
Alprazolam	309	<b>205</b>	45	42
		274		24
D5-Alprazolam	314	210	64	40
		<b>286</b>		28
Flunitrazepam	314	183	45	50
		239		34
D7-Zolpidem	315	<b>242</b>	58	40
		270		25
Clonazepam	316	214	52	40
		<b>270</b>		24
Bromazepam	316	182	24	32
		209		26
Diclazepam	319	154	40	30
		227		30
D4-Clonazepam	320	218	52	36
		<b>274</b>		22
Lorazepam	321	229	40	30
		<b>275</b>		20
8-Aminoclonazepam	324	146	54	30
		220	79	38
α-OH-alprazolam	325	216	42	38
		<b>297</b>		26
Midazolam	326	223	4	36
		291		26
Flualprazolam	327	223	62	40
		292		26
D5-α-OH-alprazolam	330	221	62	38
		<b>302</b>		24
OH-ethyl-flurazepam	333	109	52	26
		211		36
Flubromazepam	333	184	12	28
		226		32
Lormetazepam	335	177	39	42
		227		36

**Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs**

Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	13 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*



**Wisconsin Department of Justice  
Toxicology Unit Procedure Manual  
TXPM 6.5 Benzodiazepines/Z-drugs**

<b>Triazolam</b>	343	239	45	40
		308		26
<b>Etizolam</b>	343	138	6	36
		<b>314</b>		24
<b>Phenazepam</b>	348	179	45	50
		206		34
<b>Adinazolam</b>	352	205	46	50
		295		22
<b>Bromazolam</b>	353	325	40	24
		205		44
<b>Pyrazolam</b>	353	167	45	34
		206		30
<b>Clonazolam</b>	354	280	10	34
		326		26
<b>α-OH-etizolam</b>	359	315	2	20
		282		24
<b>α-OH-triazolam</b>	359	277	12	34
		331		28
<b>α-OH-clonazolam</b>	370	296	36	36
		342		26
<b>Flubromazolam</b>	371	223	12	44
		292		26
<b>α-OH-flubromazolam</b>	387	223	50	48
		359		28
<b>Flurazepam</b>	388	288	40	24
		315		24
<b>Zopiclone</b>	389	217	22	36
		245		14

Notes: Precursor ions, product ions, cone voltages and collision energies may vary

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	14 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*



# Wisconsin Department of Justice Toxicology Unit Procedure Manual TXPM 6.5 Benzodiazepines/Z-drugs

Table 9 - Quantitative analytes calibration models and weighting

Drug or Metabolite	Calibration	Weighting	Origin
<b>7-aminoclonazepam</b>	quadratic	1/x	exclude
<b>Alprazolam</b>	linear	1/x <sup>2</sup>	exclude
<b>Clonazepam</b>	quadratic	1/x	exclude
<b>Diazepam</b>	quadratic	1/x	exclude
<b>Etizolam</b>	linear	1/x <sup>2</sup>	exclude
<b>Lorazepam</b>	quadratic	1/x	exclude
<b>Nordiazepam</b>	linear	1/x	exclude
<b>Oxazepam</b>	quadratic	1/x <sup>2</sup>	exclude
<b>Temazepam</b>	quadratic	1/x <sup>2</sup>	exclude
<b>Zolpidem</b>	linear	1/x	exclude
<b>a-OH-alprazolam</b>	linear	1/x	exclude

#### 6.5.4.11 Limitations of LC/MS/MS Analysis

- No known matrix interferences were detected during validation.
- No drug interferences were detected during validation. Post-validation discoveries:
  - 4-chloro-deschloroalprazolam can be a potential interference of alprazolam.
  - Use of the precursor ion 309 and daughter 165 for alprazolam will help differentiate these analytes.
  - Further work will be conducted to add 4-chloro-deschloroalprazolam to the scope of testing, if needed.
- Quantitative results may be reported for 7-aminoclonazepam, alprazolam, clonazepam, diazepam, lorazepam, nordiazepam, oxazepam, temazepam, zolpidem, and alpha-hydroxy alprazolam at a dilution.
- Quantitative results may not be reported for etizolam in samples that are analyzed at a dilution (i.e., less than 500 µL sample volume).
- Zopiclone and eszopiclone are stereoisomers and are not differentiated via LC/MS/MS.

#### 6.5.4.12 LC/MS/MS Stability

- Validation showed that the analytes listed in table 8 are stable once extracted and stored in the cooled autosampler for **two** days in blood.

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	15 of 17



# Wisconsin Department of Justice

## Toxicology Unit Procedure Manual

### TXPM 6.5 Benzodiazepines/Z-drugs

- Bromazolam, a-OH-clonazolam, diclazepam, flualprazolam, flurazepam, phenazepam, and adinazolam are stable for **one** day in blood.
- Validation showed that the analytes listed in table 8 are stable once extracted and stored in the cooled autosampler for **two** days in urine.
  - 7-aminoclonazepam, 8-aminoclonazolam, pyrazolam, zopiclone, and 7-aminoflunitrazepam are stable for **one** day in urine.

#### 6.5.5 Results, Conclusions, and Data Analysis

These procedures are a means of extracting the substance of interest from a sample matrix in preparation for confirmation and/or quantitation by GC/MS or LC/MS/MS. See [TXPM 4.5 Gas Chromatographic and Mass Spectral Quality Control](#) or [TXPM 4.9 Liquid Chromatographic and Tandem Mass Spectral Quality Control](#) for information about the types of results and conclusions that can be made following these analyses.

- If zopiclone is confirmed, it will be reported as zopiclone/eszopiclone.
- Blood and urine samples may be analyzed by LC/MS/MS after their stability day and must be evaluated qualitatively.
  - If the analyte is acceptable, then report qualitatively. If the analyte is unacceptable, then report INC. Remediation includes reanalysis.

#### 6.5.6 Reporting

Reporting will be done after GC/MS or LC/MS/MS analysis according to the guidelines in [TXPM 3.6 Gas Chromatograph with Mass Spectrometry](#) or [TXPM 3.9 Liquid Chromatography with Tandem Mass Spectrometry](#) and TXPM 4 – Quality Guidelines.

#### 6.5.7 References

- 1) Oehldrich J & Collins CD, unpublished work, Wisconsin State Crime Laboratories
- 2) Jochemsen R & Breimer D, Journal of Chromatography, 223, p.438, 1981
- 3) Benzodiazepines and GHB Detection and Pharmacology ed. Salvatore J. Salamone: p.38 Flunitrazepam extraction pH 9-9.5 borate buffer with Diethyl Ether:Methylene Chloride (2:1) or Diethyl Ether:Chloroform (4:1); p. 58
- 4) United Chemical Technologies, Inc., Clonazepam & 7-Aminoclonazepam in Urine for GC/MS Confirmations using 200 mg Clean Screen Extraction Column
- 5) United Chemical Technologies, Inc., Benzodiazepines in Blood, plasma/serum, and tissue for LC/MS/MS Confirmations using 200 mg Clean Screen column
- 6) Colorado Bureau of Investigation TOX 10-19 Benzodiazepines and Z-Drugs quantitation by LC/MS/MS (Document #15967, Revision 8, Issue date 12/15/2022)
- 7) Milwaukee County Medical Examiner's Benzodiazepine Confirmation/Quantitation method
- 8) Wisconsin State Lab of Hygiene: Environmental Health Division, Forensic Toxicology Program, Benzodiazepines and Z-drugs Quantitation and Confirmation (revision 2.0, effective 7/12/2021)
- 9) Virginia Department of Forensic Science Benzodiazepines, Zolpidem, Zopiclone, and Zaleplon Quantification and Confirmation by LC/MS/MS (Qualtrax ID 2816, Revision 19, Issue date 6/22/2021)
- 10) Waters Application Note: Quantitative Analysis of 21 Benzodiazepine Drugs, Zolpidem and Zopiclone in Serum Using UPLC-MS/MS
- 11) Procedure For the Selection and Validation of a Calibration Model I-Description and Application, *Journal of Analytical Toxicology*, 2017;41:261–268 doi: 10.1093/jat/bkx001

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	16 of 17

This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.





# Wisconsin Department of Justice Toxicology Unit Procedure Manual TXPM 6.5 Benzodiazepines/Z-drugs

---

## 6.5.8 Safety

Universal safety precautions should be taken for biological specimens and organic waste. No novel safety concerns.

## 6.5.9 Method Validation

The GC/MS method was established prior to the Wisconsin Crime Laboratories ISO 17025 accreditation in 2012. Use of this method throughout the years has demonstrated fitness for its intended purpose.

LC/MS/MS method was approved. Files are stored in the Unit Quality Records. LC/MS/MS method validated to ANSI/ASB Standard 036, First edition, 2019 Standard Practices for Method Validation in Forensic Toxicology.

Toxicology Unit Procedure Manual — TXPM 6.5 Benzodiazepines/Z-drugs				
Document ID	Revision	Effective Date	Status	Page Number
1097	9	4/21/2023 7:14:03 AM	Published	17 of 17

*This controlled document was approved and issued in accordance with the document management provisions established in the Wisconsin Crime Laboratory Quality System. If viewed outside of Qualtrax, this is an uncontrolled copy and may no longer be the current revision.*